




WPI

ECG Classification

Yudong Yu





- **Background**
- **Methodology**
- **Analysis**





Background

- The electrocardiogram (ECG) is a diagnostic tool that is routinely used to assess the electrical and muscular functions of the heart.



Training set:
69 / 31


Test set:
64 / 36



Methodology

VGG16 + SVM


Fine-tuning VGG16



1. Feature Extractor


VGG16

a. VGG16(include_top=False) + Fully-connected layer(128)



1. Feature Extractor


- a. VGG16(include_top=False) + Fully-connected layer(128)
- b. **Transform** time series signal to **image** (RP + GASF + MTF)





1. Feature Extractor

- a. VGG16(include_top=False) + Fully-connected layer(128)
- b. **Transform** time series signal to **image** (RP + GASF + MTF)



Comparison on ECG200







2. Dimensionality Reduction



- PCA:
- n components: 20
- t-SNE: 2




3. Train SVM



- Grid search
- C: 3, Gamma: 0.01
- Best acc: 0.84
- Training ...
- Test acc: 0.80

Fine-tuning VGG16



VGG16(include_top=False)+Dense(128)+Dense(2)

L2 regularization, 0.01

BatchNormalization()

Dropout(0.5) for dense_1


Learning Rate Scheduler()

0.0002 \rightarrow 0.5*LR \rightarrow 0.2*LR



Freeze parameters in first three blocks

Training set(Augmented): 800(100*2³)



Test set: 100



Analysis(1)



FN: 92



FP: 55, 75

* these numbers are the indexes of samples in test set (100 samples)

Analysis(1)

Predictions	model 1	model 2	model 3	model 4	model 5	model 6
FN	[0, 15, 54]	[54]	[]	[]	[]	[]
FP	[9, 41]	[4, 9, 19, 22, 37]	[9, 19, 57]	[9]	[9, 19]	[9, 19]
False in total	5	6	3	1	2	2
Test Acc.	0.92	0.91	0.93	0.95	0.93	0.94

* these number are the indexes of samples in test set(60 samples)

Analysis(1)

Predictions	model 1	model 2	model 3	model 4	model 5	model 6
FN	[0, 15, 54]	[54]	[]	[]	[]	[]
FP	[9, 41]	[4, 9, 19, 22, 37]	[9, 19, 57]	[9]	[9, 19]	[9, 19]
False in total	5	6	3	1	2	2
Test Acc.	0.92	0.91	0.93	0.95	0.93	0.94

Analysis(2)

Table 1. Accuracy comparison with traditional classification algorithm(The method/s with the highest accuracy in each database are shown in bold)

(a)						
Data	C4.5	C4.5(S)	1NN	1NN(S)	NaB	NaB(S)
Adiac	53.19	49.36	59.34	56.27	56.52	57.54
Beef	56.67	40	60	53.33	50	60
Chlorine	64.3	56.82	68.52	58.59	34.61	45.52
Coffee	57.14	92.86	75	100	67.86	92.86
Diatom	71.24	67.65	93.46	94.44	87.91	78.76
ECG200	72	79	89	78	77	80

2016

Table 2 Error rates of the methods for a subset of UCR dataset

Dataset	DTW	ST	TSBF	HOG1D	CovNN	CovSVM
50 Words	0.242	0.281	0.209	0.402	0.222	0.200
Adiac	0.391	0.435	0.245	0.320	0.217	0.164
Beef	0.467	0.167	0.287	0.367	0.100	0.067
CBF	0.004	0.003	0.009	0.000	0.000	0.000
ChlorineCon	0.350	0.300	0.336	0.307	0.294	0.255
CinCECGT	0.070	0.154	0.262	0.249	0.003	0.000
Coffee	0.179	0.000	0.004	0.000	0.000	0.000
CricketX	0.236	0.218	0.278	0.195	0.244	0.236
CricketY	0.197	0.236	0.259	0.205	0.210	0.251
CricketZ	0.180	0.228	0.263	0.185	0.239	0.215
DiatomSizeR	0.065	0.124	0.126	0.016	0.052	0.043
ECG200	0.310		0.145	0.060	0.080	0.070

June 2018





Table 5. Overall classification accuracies of *k*-NN alignment with six scenarios: NONE, DTW, CDTW, SAGA, PTW and CTW.

Dataset	NONE	DTW	CDTW	SAGA	PTW	CTW
Synthetic Control	88.0	99.3	98.7	90.3	94.3	98.7
Trace	76.0	100.0	99.0	99.0	99.0	100.0
Sony AIBO Robot Surface II	85.9	83.1	85.9	86.5	85.0	84.3
Sony AIBO Robot Surface	69.6	72.5	69.6	73.5	74.7	73.0
Symbols	89.9	95.0	93.8	94.9	95.3	90.1
Two Lead ECG	74.7	90.4	86.8	86.7	98.4	90.3
Olive Oil	86.7	86.7	83.3	83.3	86.7	80.0
Mote Strain	87.9	83.5	87.9	89.9	86.6	84.0
Lighting 7	57.5	72.6	71.2	79.5	64.4	68.5
Lighting 2	75.4	86.9	86.9	85.2	80.3	85.2
Italy Power Demand	95.5	95.0	95.5	96.4	94.5	95.0
Gun Point	91.3	90.7	91.3	98.7	98.7	88.0
Face Four	78.4	83.0	88.6	77.3	87.5	85.2
ECG Five Days	79.7	76.8	79.7	96.9	90.5	75.3
ECG 200	88.0	77.0	88.0	87.0	85.0	79.0

2015

Analysis(3)





RP + GASF + MTF



Processing
1D time series signal
with **transformation
methods** or not



With a simple NN
dense(50)+dense(2)

ECG5000



Question?

ECG200



ECG5000

